



Bellwether Porcelain Crucible Lids

Description • Lids only

- To assure quality, the porcelain is oven heated to 1100 DEG C, cooled to 400 DEG C in oven, removed from oven and cooled to 20 DEG C three times

Image(s) are representative of the product group and not necessarily the individual product.

Order	29-Apr-2024 04:53: AM PDT	PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.					
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
980-40969-EA	36 mm I.D. Porcelain Crucible Lid - 1 ea Lids Use with: 980-40972 	Sign In	Vee Gee Scientific	5011L			EA
980-40969-EA	36 mm I.D. Porcelain Crucible Lid - 1 ea Lids Use with: 980-40972 		Vee Gee Scientific	5011L			CS(10)
980-41014-EA	71 mm I.D. Porcelain Crucible Lid - 1 ea Lids Use with: 980-41017 		Vee Gee Scientific	5015L			EA
980-41014-EA	71 mm I.D. Porcelain Crucible Lid - 1 ea Lids Use with: 980-41017 		Vee Gee Scientific	5015L			CS(10)