

Beakers for Bransonic* Ultrasonic Cleaning Baths



Description • Designed to maximize the effectiveness of Bransonic* baths

- It is important that parts do not rest on the bottom of the bath, as it can inhibit the ultrasonic activity and possibly damage the unit
- Indirect cleaning method is ideal for sample preparation
- Suspends parts using beakers with a beaker positioning cover†
- Ideal for cleaning small parts
- Discarded soil cleans easily from a beaker or solid tray

†Available separately, see "Beaker Positioning Covers for Bransonic* Ultrasonic Baths"

Image(s) are representative of the product group and not necessarily the individual product.

| Order | 28-Apr-2024 13:49: PM PDT | | PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE. | | | | |
|--------------|---|-------------------------|--|-------------|------|------------|------|
| Catalog | Item | Availability | Manufacturer | Mfg. Part # | NDC# | Your Price | Unit |
| 382-10245-EA | 250 mL Glass Beaker for use with Ultrasonic Cleaner - 1 ea Material: Glass Capacity (mL): 250 | Sign In | Branson Ultrasonics Corporation | 000-140-001 | | | EA |