

Burkle TargetDispo LaboPlast and SteriPlast Disposable Samplers



Description TargetDispo is used to take a target point sample at a position in the sample material that has previously been precisely determined. The sampler pierces into the powder or granulate while closed. When it reaches the sample removal point it is opened and then closed again. In this way a precise target-point sample is obtained. The sample can simply be emptied through the open handle.

- Easy emptying through the opening in the handle
- HDPE
- Cleanroom manufactured
- Optionally sterilized by gamma rays

Image(s) are representative of the product group and not necessarily the individual product.

No options of this product are available.