



WWW.SPECTRUMRX.COM



Description • Reduces the possibility of accidentally spilling liquid waste

- Simple four step process takes less than two minutes
- Disposal bags and LabSorb granules can be used in conjunction with most 500 mL, 600 mL, 1 L and 2 L beakers
- LabSorb granules can rapidly absorb aqueous or mixed aqueous and organic solvent waste
- Convenient system includes extra durable disposal bags, super absorbent LabSorb granules and a convenient graduated scoop
- •When spills do occur, use LabSorb for safe liquid containment during cleanup†
- Use at locations connected to septic systems to avoid pouring waste down the drain
- Put Lab-Sorb at the bottom of your biohazard bin
- Convert liquids to solids in seconds
- One jar of LabSorb can absorb up to 50 liters of liquid
- Retains liquid during autoclaving

†Always follow proper safety protocols

Image(s) are representative of the product group and not necessarily the individual product.

Order	29-Apr-2024 01:31: AM PDT	PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.					
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
110-51010- EA	Diversified Biotech Lab-Sorb Granules, 16 oz Jar Only - 1 ea Lab-Sorb Only Bag Size: n/a Bag Dimensions (in.): n/a	Sign In	Diversified Biotech, Inc.	SORB- 1000			EA
110-51010- EA	Diversified Biotech Lab-Sorb Granules, 16 oz Jar Only - 1 ea Lab-Sorb Only Bag Size: n/a Bag Dimensions (in.): n/a		Diversified Biotech, Inc.	SORB- 1000			CS(5)

