

EM Science Histology Alcohols and Solvents

Description • EM Science Histology Alcohols and Solvents. Used to hydrate or dehydrate tissue samples in tissue processing and staining protocols

- Additional shipping charges will be added for hazardous items



- Harleco Diluted Alcohols
- For histology and cytology, these alcohols are available at 100%, 95%, 80% and 70% concentrations in easy-to-pour poly bottles
- Neo-Clear Xylene Substitute
- Aromatic solvent used as a medium for clearing and deparaffinizing steps in staining and tissue processing
- HARLECO Acid Alcohol
- Decolorizer for differentiating nuclear chromatin in regressive staining methods

Image(s) are representative of the product group and not necessarily the individual product.

No options of this product are available.

Shipping & Storage

Shipping Information

DOT: UN1986, Alcohols, Flammable, Toxic, N.O.S. (Cont. Methanol, Ethanol, Propanol), 3, PG II