

SP Bel-Art Colony Counter

Description • A cost effective and accurate way to count colonies of bacteria or other microorganisms grown on an agar plate

- This simple to operate device marks, counts, and gives an audible and visual alert each time a count is registered
- The lightweight, pen-style design reduces the likelihood of error, combining all aspects of counting into one motion
- F37862-0002 is CE marked



- Fast, accurate and economical sample counting and sample count verification
- Stores up to 35 sample counts
- Verify sample counts with up/down counting
- Displays sum total of all colonies counted
- Backlit display for easy viewing in low light conditions
- Audible alert can be turned off if desired
- Counting function can be locked to prevent error during breaks between counting
- Long life, red and black felt-tip markers (included) work on plastic or glass and are oil and water resistant; ink is removable with ethyl alcohol
- Includes 3V lithium battery and protective storage case

Image(s) are representative of the product group and not necessarily the individual product.

Order	28-Mar-2024 14:54: PM PDT	PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.					
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
144-78514-EA	SP Bel-Art Colony Counter - 1 ea	Sign In	SP Bel-Art	F37862-0002			EA
144-78514-EA	SP Bel-Art Colony Counter - 1 ea		SP Bel-Art	F37862-0002			CS(4)