

## C & G VOA Vials with HCl Preservative

**Description** • VOA Vials are most commonly used for Volatile Organic Compound Analysis

- These vials are Type 1 Borosilicate glass. Vials come with a 2pc. 0.125 in. PTFE/Silicon septum press fitted into 24-414 white hole closure

Image(s) are representative of the product group and not necessarily the individual product.

Order	05-May-2021 21:42: PM PDT	PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.					
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
980-68100-CS	60ml Amber, VOA Vial w/ 6mg of Ammonium Chloride, Pre-cleaned, Certified, 72/case - 72 ea	<a href="#">Sign In</a>	Spectrum Chemical	LY84-E110-C61A			CS
980-68103-CS	40ml Clear, VOA Vial w/ 3 mg of Sodium Thiosulfate, Precleaned, Certified - 72 ea		Spectrum Chemical	LV74-B310-C31A			CS

05-May-2021 21:42: PM PDT