

## C & G VOA Vials with HCl Preservative



**Description** • VOA Vials are most commonly used for Volatile Organic Compound Analysis

- These vials are Type 1 Borosilicate glass. Vials come with a 2pc. 0.125 in. PTFE/Silicon septum press fitted into 24-414 white hole closure

Image(s) are representative of the product group and not necessarily the individual product.

Order	28-Mar-2024 10:19: AM PDT		PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.				
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
980-68103-CS	40ml Clear, VOA Vial w/ 3 mg of Sodium Thiosulfate, Precleaned, Certified - 72 ea	<a href="#">Sign In</a>		LV74-B310-C31A			CS