



## C & G VOA Vials

**Description** • VOA Vials are most commonly used for Volatile Organic Compound Analysis

- These vials are Type 1 Borosilicate glass available in Clear and Amber
- Vials come with a choice of 1pc. 24-414 white cap and sonically bonded septum, or 2pc. 0.125 in. PTFE/Silicone septum press fitted into 24-414 white hole closure

Image(s) are representative of the product group and not necessarily the individual product.

Order	25-Apr-2024 00:59: AM PDT		PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.				
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
980-68174-CS	40ml Flint Glass Vial with 24-414 White 2 piece PTFE Septa cap level B-1/COV - 72 ea	<a href="#">Sign In</a>		LV73-A010-A01B			CS