

## Veltek DECON-AHOL WFI\* Formula Sterile 60% Isopropyl Alcohol (IPA)

[SDS](#)



**Description** • VAI's DECON-AHOL\* Products are formulated with USP Isopropyl Alcohol (IPA)

- DECON-AHOL 60% USP IPA with USP Water for Injection (WFI) has been designed specifically for pharmaceutical, biotechnology, healthcare, and medical device cleaning rotations for use in both aseptic and non-aseptic environments
- DECON-AHOL WFI Formula can be used for the decontamination of numerous surfaces within any cleanroom operation
- Each lot of DECON-AHOL is sterility tested according to current USP Compendium, is completely traceable, and has been completely validated for sterility and shelf life
- One gallon containers are an excellent choice for GMP operations that require a large volume of solution

- Raw components are air washed with 0.2 micron filtered air
- Filled in an ISO 5 (Grade A/B, Former Class 100)
- Filtered at 0.2 microns Gamma irradiated at a  $10^{-6}$  SAL
- Lot sterility tested according to current USP compendium
- Completely lot traceable from start to finish
- Completely validated for sterility and shelf life

Image(s) are representative of the product group and not necessarily the individual product.

Order	29-Apr-2024 12:47: PM PDT		PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.				
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
120-97735-CS	Veltek Associates, Inc. Veltek DECON-AHOL WFI Formula Sterile 60 Percent Isopropyl Alcohol (IPA), 1 gallon, Sterile - 4 ea	<a href="#">Sign In</a>	Veltek Associates, Inc.				CS

### Shipping & Storage

#### Shipping Information

DOT: UN1219, Isopropanol, 3, PG II