






Veltek DECON-AHOL WFI* Formula Sterile 91% Isopropyl Alcohol (IPA)

SDS



- DECON-AHOL 91% USP IPA with USP Water for Injection (WFI) has been designed specifically for pharmaceutical, biotechnology, healthcare, and medical device cleaning rotations for use in both aseptic and non-aseptic environments
- DECON-AHOL WFI Formula can be used for the decontamination of numerous surfaces within any cleanroom operation
- Each lot of DECON-AHOL is sterility tested according to current USP Compendium, is completely traceable, and has been completely validated for sterility and shelf life
- One gallon containers are an excellent choice for GMP operations that require a large volume of solution
- Raw components are air washed with 0.2 micron filtered air
- Filled in an ISO 5 (Grade A/B, Former Class 100)
- Filtered at 0.2 microns Gamma irradiated at a 10⁻⁶ SAL
- Lot sterility tested according to current USP compendium
- Completely lot traceable from start to finish
- Completely validated for sterility and shelf life

Image(s) are representative of the product group and not necessarily the individual product.

| Order | 28-Apr-2024 16:03: PM PDT | PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE. | | | | | |
|------------------|--|--|----------------------------|----------------|------|---------------|------|
| Catalog | Item | Availability | Manufacturer | Mfg. Part # | NDC# | Your Price | Unit |
| 120- 97744-CS | Veltek Associates, Inc. Veltek DECON-AHOL WFI Formula Sterile 91 Percent Isopropyl Alcohol (IPA), 1 gallon, Sterile - 4 ea | Sign In | Veltek Associates, Inc. | | | | CS |

Shipping & Storage

Shipping Information

DOT: UN1219, Isopropanol, 3, PG II