

Sartorius Stedim BACTair Empty Plates

Description • BACTair empty plates for microbiological air monitoring with the AirPort MD8 air sampler (impaction method) **Media not included**

- The system allows impaction onto the culture media plate, since the plate itself functions as collection head
- Individually, sterile packaged
- Integrated disposable sieve
- Samples 1 m³ in 8 minutes
- Optimized geometry; the number of impaction-holes (400) and their accurate positioning are key to reliable results



Specifications:

- Dimensions: 11.6 x 2.4 cm
- Size of impaction holes: 0.47 mm

Image(s) are representative of the product group and not necessarily the individual product.

Order	29-Apr-2024 13:58: PM PDT	PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.					
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
946-10011-PK	BACTair Empty Plates without media, sterile - 50 ea	Sign In	Sartorius	14301-110-----K			PK