spectrum®

RE: Residual Solvents - Finasteride, USP (Cat# F1300)

To Whom It May Concern:

Thank you for your interest in Spectrum high quality chemicals.

The above mentioned catalog item complies with the requirements of USP <467> Residual Solvents and ICH Q3C Residual Solvents Guideline. The following solvents are likely to be present:

Residual Solvent	Class	Expected Concentration
Methanol	2	< 3000 ppm
Methylene Chloride	2	< 600 ppm
Toluene	2	< 890 ppm
Ethanol	3	< 5000 ppm
Isopropyl Acetate	3	< 5000 ppm

No Class 1 or other solvents are used or produced in the manufacturing process of this material.

For lot-specific information, please refer to the respective Certificate of Analysis.

If you have any further questions, please contact us by telephone at 1(800)772-8786 Option 2, or by email at <u>TechServices@spectrumchemical.com</u>.

Sincerely,

Technical Services Spectrum Chemical Mfg. Corp.

This document has been produced electronically and is valid without a signature.

CORPORATE HEADQUARTERS

769 Jersey Avenue New Brunswick, NJ 08901 PHONE 732.214.1300 FAX 732.246.7132

WEST COAST FACILITY

14422 South San Pedro Street Gardena, California 90248 PHONE 310.516.8000 FAX 310.516.9843

AN ISO 9001:2015 REGISTERED COMPANY