

Chex-All II Instant Sealing Pouches

Description • For use in steam autoclaves or ethylene oxide sterilizers

- Provides a barrier to maintain sterility after the sterilization process
- Pouches have indicator imprints on both the inside and outside
- Monitors the effective penetration of steam or ethylene oxide gas on the inside of the package
- These self-sealing pouches do not require heat sealers
- Provides a seal as effective in maintaining sterility as any heat seal
- In steam autoclaves, the yellow indicators turn dark
- In ethylene oxide sterilizers, the red indicators change to green
- $\frac{3}{4}$ in. wide adhesive line guarantees easy and reliable sealing
- Triple sealed paper/plastic film provides excellent sterility protection
- Hermetic seal will not open in steam autoclaves and is not affected by ethylene oxide gas
- Chevron with corner seals assures dust-free, sterile presentation
- Transparent side of the pouches facilitates visual inspection of sterilizable materials and instruments
- Premium quality 2.5 mil. plastic and 73 - 74 g medical grade paper
- Resistant to tearing, bursting, fracturing, and bacterial penetration
- Provides a high degree of protection and allows for a reliable transfer to the sterile field



Image(s) are representative of the product group and not necessarily the individual product.

Order	28-Apr-2024 08:25: AM PDT	PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.					
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
311-10143-CS	Chex-All II Instant Sealing Pouches, 5 x 10 in. (12.7 x 25.4 cm) - 500 ea Dimensions W x L: 5 x 10 	Sign In		02401000			CS