

Sterile Filtration Cups Bubble Point Test Apparatus



- Description** • Ensure product integrity of filter units or bottletop filters
- Tests filter unit membranes without removing the membrane from the filter unit
 - Cast aluminum
 - Tests the bubble point of 50, 75, and 90 mm filtration products
 - Calibrate gauge and pressure source not included
 - Accepts 0.25 in. (6 mm) I.D. tubing
 - Includes adapters and O-rings

Image(s) are representative of the product group and not necessarily the individual product.

Order	29-Apr-2024 05:42: AM PDT		PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.				
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
147-20194-EA	Thermo Scientific Nalgene Bubble Point Test Apparatus For all 50, 75 and 90 mm MF75 filter units and bottle top filters - 1 ea	Sign In	Thermo Scientific Nalgene	DS0405-0050			EA