







Labconco 6-Port Manifold for FreeZone Bulk Tray Dryers and Stoppering Tray Dryers

Description • Six valves provide the flexibility to connect flasks and other freeze dry glassware to the Freeze Dry System with Stoppering or Bulk Tray Dryer

- Valves accommodate either $\frac{1}{1/2}$ in. or $\frac{3}{4}$ in. flask adapters
- Mounts to support stand included with the Tray Dryer
- May be field installed before or after the Tray Dryer is installed on the Freeze Dry System
- Included with Bulk Tray Dryer
- For use with the FreeZone Triad (Cat. No. 416-10001), FreeZone Stoppering Tray Dryer (Cat. No. 346-27512) and FreeZone Bulk Tray Dryer (Cat. No. 496-29438 and 496-29441).

Image(s) are representative of the product group and not necessarily the individual product.

Order	30-Apr-2024 00:29: AM PDT	PRICES ARI	E SUBJECT TO CI	HANGE WI	THOUT P	RIOR NO	TICE.
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
496- 28058-EA	6-Port Manifold for FreeZone Bulk Tray Dryers and Stoppering Tray Dryers - 1 ea Dimensions W x D x H in. (cm): 30.0 x 16.5 x 7.0 (76.2 x 41.9 x 17.8) Shipping Weight lbs (kg): 19 (8.6) 	Sign In	Labconco	7726500			EA