

## Topi-CLICK PERL Vaginal Dosing System



**Description** • Vaginal dosing system with comfortable and accurate applicator

- The dispenser is a Topi-Click 35 (up to 37 mL capacity) with a hinged shell cap and Snap-on Applicator Port
- The Snap-on Applicator is made with FDA approved clear, flexible soft medical grade polymer
- Provides precise and accurate dosing, 0.25 mL per click
- Produces significantly less residual than the traditional tube and applicator
- Easy-to-use and easy-to-clean

Image(s) are representative of the product group and not necessarily the individual product.

| Order        | 19-Apr-2024 17:48: PM PDT                                                        | PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE. |              |                  |      |            |        |
|--------------|----------------------------------------------------------------------------------|----------------------------------------------------|--------------|------------------|------|------------|--------|
| Catalog      | Item                                                                             | Availability                                       | Manufacturer | Mfg. Part #      | NDC# | Your Price | Unit   |
| 554-10138-PK | Topi-CLICK PERL 35 mL Vaginal Dosing System with 4 mL Vaginal Applicator - 25 ea | <a href="#">Sign In</a>                            |              | TCPERL01 - WHITE |      |            | PK     |
| 554-10138-PK | Topi-CLICK PERL 35 mL Vaginal Dosing System with 4 mL Vaginal Applicator - 25 ea |                                                    |              | TCPERL01 - WHITE |      |            | CS(10) |
| 554-10141-PK | Topi-CLICK PERL 35 mL Vaginal Applicator Dose Loader - 25 ea                     |                                                    |              | TCPERL02         |      |            | PK     |
| 554-10141-PK | Topi-CLICK PERL 35 mL Vaginal Applicator Dose Loader - 25 ea                     |                                                    |              | TCPERL02         |      |            | CS(10) |
| 554-10144-PK | Topi-CLICK PERL 4 mL Vaginal Dose Applicator - 25 ea                             |                                                    |              | TCPERL03         |      |            | PK     |
| 554-10144-PK | Topi-CLICK PERL 4 mL Vaginal Dose Applicator - 25 ea                             |                                                    |              | TCPERL03         |      |            | CS(10) |