



## Biospec Chrome-Steel Beads

**Description** • Chrome steel has a density of 7.9g/cc

- Used mostly for grinding leaves
- Beads are sold in one pound bottles
- With chrome steel beads, common polypropylene microvials will break if 'dry' grinding or cryo-grinding is attempted; stainless steel microvials with silicone rubber caps or XXTuff poly vials are recommended

Image(s) are representative of the product group and not necessarily the individual product.

Order	29-Apr-2024 03:56: AM PDT		PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.				
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
140-20059-EA	Chrome-Steel Beads, 3.2 mm diameter - 1 ea Diameter (mm): 3.2 	<a href="#">Sign In</a>	Biospec Products	11079132c-B			EA
140-20077-EA	Chrome-Steel Beads, 6.35 mm diameter - 1 ea Diameter (mm): 6.35 		Biospec Products	11079635c-B			EA