

## Labconco Attachment Port Lid



**Description** • Includes  $\frac{3}{4}$  in. thick, clear acrylic cover with 3 in. diameter opening and gasket and acrylic plus with gasket (for periods of non-use)  
 • Lid replaces the lid to the built-in drying chamber on FreeZone 4.5 Liter Freeze Dry Systems to permit use of the Heated Drying Chamber (Cat. No. 336-21913), 4-Port Manifold (Cat. No. 346-28361), 20-Port Manifold (Cat. No. 346-17191), 48-Port Manifold (Cat. No. 346-10084), or Clear Chamber (Cat. No. 336-23014)

Image(s) are representative of the product group and not necessarily the individual product.

Order	29-Apr-2024 12:01: PM PDT	PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.					
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
346-17719-EA	10.38 in. Attachment Port Lid for 4.5L Freeze Dryers - 1 ea	<a href="#">Sign In</a>	Labconco	7762800			EA