

714083

Description • The job of collecting water samples is easier with the Whirl-Pak Stand-Up Thio-Bag

- A special gusset on the bottom of the bag flattens out when water is put in the bag, allowing it to stand up completely by itself
- No rack or other holder is needed to keep the bag upright
- Forget washing and sterilizing glass bottles and disposing of expensive vials
- Sterile and unbreakable, these bags will save time and money whether working in the lab or in the field
- Each bag contains 10 mg of active sodium thiosulfate to neutralize chlorine in the sample
- Capacity is measured in fluid ounces when the bag is closed and tab is folded over three times
- Volume and dimensions are approximate
- Bags should not be used at temperatures above 180°F (82°C)
- Bags can be frozen to any temperature, even liquid nitrogen temperatures down to -346°F (-210°C) - careful handling required after freezing
- All bags are sterilized after manufacturing



Image(s) are representative of the product group and not necessarily the individual product.

No options of this product are available.