



Oakton Acorn All-in-One pH/ATC Electrodes

Description • All-in-one pH/ATC Probes are sealed, have an epoxy bodd and a 3-ft cable with Sub-mini/BNC connector

- 12 mm diameter and 115 mm length
- Use with Oakton pH meters accepting 2.5 mm audio plug for temperature
- Available with single or double junction

Image(s) are representative of the product group and not necessarily the individual product.

Order	27-Apr-2024 03:05: AM PDT		PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.				
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
244-11179-EA	Oakton All-in-One pH/ATC Electrode, Epoxy, Sealed, Double-Junction, 12 mm OD; 3 ft Sub-Mini/BNC - 1 ea Junction: Double 	Sign In	Antylia Scientific	WD-35811-72			EA
244-27370-EA	Oakton All-in-One pH/ATC Electrode, Epoxy, Sealed, Single-Junction, 12 mm OD; 3 ft Sub-Mini/BNC - 1 ea Junction: Single 		Antylia Scientific	WD-35811-71			EA