





## **3M ATTEST 48-Hour Biological Indicators** for Steam

**SDS** 

**Description** • Biological Indicator for steam 270°F/132°C vacuum assisted or 250°F/121°C gravity displacement sterilization cycles, 48 hour readout, Brown cap

- 3M Attest Biological indicators that monitor the effectiveness of the steam sterilization process and provide results within 48 hours
- The self-contained BI for 250°F/121°C gravity and 270°F/132°C vacuum-assisted steam cycles consist of a Geobacillus stearothermophilusspore strip; sealed glass ampule with growth medium and a bromocresol purple pH indicator system; brown color-coded cap with holes for sterilant penetration and a filter that acts as a barrier, and a process indicator on the label that changes from rose to brown when processed
- After sterilization, the vial is crushed to join the growth media with the processed spore strip
- The BI is incubated for 48 hours for a visual change readout. A color change to yellow indicates surviving spores and a positive result.
- Biological indicator vials are easy to use and interpret with visual color readout in 48 hours
- Attest Biological Indicator 1262 for 250°F/121°C Gravity and 270°F/132°C Vacuum Assisted Steam Sterilizers
- Use in conjunction with 3M Attest Biological Indicator Incubator 116
- Can eliminate recalls when used in every load
- 1262-S is designed for industrial sterilization applications (Not FDA Approved Life Science Labs Only)

†Continental US+HI Only - No Exports

Image(s) are representative of the product group and not necessarily the individual product.

Order	24-Apr-2024 22:49: PM PDT	PRICES ARE	SUBJECT TO CH	ANGE W	ITHOUT	PRIOR NO	OTICE.
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
312- 10145-CS	3M ATTEST Biological Indicators for Steam Sterilization, 48-Hour Readout, Brown Cap, 100/box, 4 boxes/case - 400 ea Units per Box: 100 Sterility: Non- Sterile 	Sign In	3M				CS
312- 10148-CS	3M ATTEST Biological Indicators for Steam Sterilization, 48-Hour Readout, Brown Cap, 25/box, 4 boxes/case - 100 ea Units per Box: 25 Sterility: Non-Sterile 		3M				CS
312- 10151-CS	3M ATTEST Biological Indicators for Steam Sterilization, Sterile, 48-Hour Readout, Brown Cap, 300/box, 2 boxes/case - 600 ea Units per Box: 300 Sterility: Sterile 		3M				CS



## Shipping & Storage

## **Storage Information**

59-86?F (15-30?C), 35-60% relative humidity