

Topi-CLICK MICRO Topical Dosing Applicators, Round Applicator

Description The Topi-CLICK MICRO is a revolutionary way to accurately and easily micro-dose topical medications for humans and pets

Accurate & Consistent Dose

- Topi-CLICK dispensers proven most accurate among others tested by ARL
- Eliminates user-variability
- Consistent automated dosing throughout treatment

Automated-Metered Dosing

• Audible & tactile and visual CLICKs indicate measure of dose at each point of the unique triad shaped base

• One CLICK = 0.05 mL

Targeted Topical Therapy

- Innovative triad shaped dispenser allows for precision handling during application
- Reduces the risk of hand transference to others with applicator
- Offers applicator selections for customized dosing

Viscosity Assurance

- Accuracy NOT affected by varying formulations
- · Accuracy NOT affected by varying bases
- Dispense per CLICK: 0.05 mL
- Delivery volume: 1 mL to 9.75 mL
- Approx. Residual: 0.15 mL
- Made of FDA grade materials
- Phthalate free
- UV block and is recyclable

Image(s) are representative of the product group and not necessarily the individual product.

| Order | 18-Apr-2024 08:22: AM PDT | PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE. | | | | | |
|------------------|--|--|--------------|----------------|------|---------------|--------|
| Catalog | Item | Availability | Manufacturer | Mfg. Part # | NDC# | Your Price | Unit |
| 554-10195- PK | Topi-CLICK Micro with Rounded Applicator, Blue, 9 mL (0.05 mL/click) - 25 ea | <u>Sign In</u> | | TC9-R- Blue | | | РК |
| 554-10195- PK | Topi-CLICK Micro with Rounded Applicator, Blue, 9 mL (0.05 mL/click) - 25 ea | | | TC9-R- Blue | | | CS(10) |

Spectrum[®]