

Topi-CLICK MICRO 9 mL Topical Dosing Applicators, Blue

Description The Topi-CLICK MICRO is a revolutionary way to accurately and easily micro-dose topical medications for humans and pets

Accurate & Consistent Dose

- Topi-CLICK dispensers proven most accurate among others tested by ARL
- Eliminates user-variability
- Consistent automated dosing throughout treatment

Automated-Metered Dosing

- Audible & tactile and visual CLICKs indicate measure of dose at each point of the unique triad shaped base
- One CLICK = 0.05 mL

Targeted Topical Therapy

- Innovative triad shaped dispenser allows for precision handling during application
- Reduces the risk of hand transference to others with applicator
- Offers applicator selections for customized dosing

Viscosity Assurance

- Accuracy NOT affected by varying formulations
- Accuracy NOT affected by varying bases

- Dispense per CLICK: 0.05 mL
- Delivery volume: 1 mL to 9.75 mL
- Approx. Residual: 0.15 mL
- Made of FDA grade materials
- Phthalate free
- UV block and is recyclable



Image(s) are representative of the product group and not necessarily the individual product.

Order	18-Apr-2024 06:29: AM PDT	PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.					
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
554-10189-PK	Topi-CLICK Micro with Soft Angled Applicator, Blue, 9 mL (0.05 mL/click) - 25 ea	Sign In		TC9-A-Blue			PK
554-10189-PK	Topi-CLICK Micro with Soft Angled Applicator, Blue, 9 mL (0.05 mL/click) - 25 ea			TC9-A-Blue			CS(10)
554-10192-PK	Topi-CLICK Micro with Pin Point Applicator, Blue, 9 mL (0.05 mL/click) - 25 ea			TC9-P-Blue			PK

Order	18-Apr-2024 06:29: AM PDT	PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.					
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
554-10192-PK	Topi-CLICK Micro with Pin Point Applicator, Blue, 9 mL (0.05 mL/click) - 25 ea			TC9-P-Blue			CS(10)
554-10195-PK	Topi-CLICK Micro with Rounded Applicator, Blue, 9 mL (0.05 mL/click) - 25 ea			TC9-R-Blue			PK
554-10195-PK	Topi-CLICK Micro with Rounded Applicator, Blue, 9 mL (0.05 mL/click) - 25 ea			TC9-R-Blue			CS(10)