

WWW.SPECTRUMRX.COM

# **Topi-CLICK MICRO 9 mL Topical Dosing Applicators, White**

**Description** The Topi-CLICK MICRO is a revolutionary way to accurately and easily micro-dose topical medications for humans and pets

## **Accurate & Consistent Dose**

- Topi-CLICK dispensers proven most accurate among others tested by ARL
- Eliminates user-variability
- Consistent automated dosing throughout treatment

### **Automated-Metered Dosing**

- Audible & tactile and visual CLICKs indicate measure of dose at each point of the unique triad shaped base
- One CLICK = 0.05 mL

# **Targeted Topical Therapy**

- Innovative triad shaped dispenser allows for precision handling during application
- Reduces the risk of hand transference to others with applicator
- Offers applicator selections for customized dosing

### **Viscosity Assurance**

- Accuracy NOT affected by varying formulations
- Accuracy NOT affected by varying bases
- Dispense per CLICK: 0.05 mL
- Delivery volume: 1 mL to 9.75 mL
- Approx. Residual: 0.15 mL
- Made of FDA grade materials
- · Phthalate free
- UV block and is recyclable

Image(s) are representative of the product group and not necessarily the individual product.

| Order            | 26-Apr-2024 10:28: AM PDT   | PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE. |              |                |      |               |        |
|------------------|---|--|--------------|----------------|------|---------------|--------|
| Catalog          | Item  | Availability                                       | Manufacturer | Mfg.<br>Part # | NDC# | Your<br>Price | Unit   |
| 196-13867-<br>PK | Topi-CLICK Micro with Soft Angled<br>Applicator, White, 9 mL (0.05 mL/click) - 25<br>ea | Sign In  |              |                |      |               | РК     |
| 196-13867-<br>PK | Topi-CLICK Micro with Soft Angled<br>Applicator, White, 9 mL (0.05 mL/click) - 25<br>ea |  |              |                |      |               | CS(10) |

