

California Pacific Labs Luer Barb HPLC Waste Line Adapters

Description • ECO Funnels can be equipped with HPLC fittings, allowing the operator to connect HPLC waste lines directly to the funnel

- Installation is free of charge. To order HPLC attached to your ECO Funnel(s), please add the desired funnels and fittings to your cart
- We can install up to 6 HPLC fittings per 4 in. funnel and 8 fittings per 8 in. funnel
- Measure the Inner Diameter (ID) of your tubing to determine the size of Luer Barb fitting you need

- These luer fittings, luer adapters and luer accessories are made to the ISO-594-1 standard, and may be connected to other luers meeting ISO standards
- Please test your application to verify the performance of the interconnection
- Part of a standardized system of small-scale tubing connectors used to make secure leak-proof unions in medical devices and laboratory instruments, these luer fitting provide a unique secure friction seal appropriate for gas and fluid connections


- The chemically resistant polypropylene Luer barb adapters have interconnected Luer fittings
- Tubing connections range from $\frac{1}{16}$ in. (1.6 mm) to $\frac{1}{4}$ in. (6.4 mm) The threaded end is connected to the ECO Funnel while the tubing connects to the Luer barb

- The barb has no parting lines on the barb to ensure a leak-free tubing union
- Female lock luers utilize lugs or threads that engage securely with special threads on the male lock luers
- Male lock luers are available in versions with integral lock rings or with rotatable features to permit tubing positioning prior to assembly

Each HPLC adapter includes a male and female side:

1. Female Luer Thread Style with $\frac{5}{16}$ in. Hex to 1/4-28 UNF Bottom Sealing Thread, Natural Polypropylene
2. Male Luer Integral Lock Ring to Barb, Natural Polypropylene, available in several sizes to fit different tubing:

&



Order	01-May-2024 03:04: AM PDT	PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.					
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
196-20279-EA	HPLC Kit, Luer Barb PP Adapter for 1/16 in. ID Tubing - 1 ea Size (in.): 1/16 	Sign In		HPLC-LB-1/16-KIT			EA

Order	01-May-2024 03:04: AM PDT	PRICES ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE.					
Catalog	Item	Availability	Manufacturer	Mfg. Part #	NDC#	Your Price	Unit
196-20280-EA	HPLC Kit, Luer Barb PP Adapter for 1/8 in. ID Tubing - 1 ea Size (in.): 1/8 			HPLC-LB-1/8-KIT			EA
196-20281-EA	HPLC Kit, Luer Barb PP Adapter for 1/4 in. ID Tubing - 1 ea Size (in.): 1/4 			HPLC-LB-1/4-KIT			EA
196-20282-EA	HPLC Kit, Luer Barb PP Adapter for 3/32 in. ID Tubing - 1 ea Size (in.): 3/32 			HPLC-LB-3/32-KIT			EA
196-20283-EA	HPLC Kit, Luer Barb PP Adapter for 5/32 in. ID Tubing - 1 ea Size (in.): 5/32 			HPLC-LB-5/32-KIT			EA
196-20284-EA	HPLC Kit, Luer Barb PP Adapter for 3/16 in. ID Tubing - 1 ea Size (in.): 3/16 			HPLC-LB-3/16-KIT			EA